

# ***Multiphoton Microscopy in the Biomedical Sciences XVIII***

**Ammasi Periasamy**

**Peter T. C. So**

**Karsten König**

**Xiaoliang S. Xie**

*Editors*

**28–30 January 2018**

**San Francisco, California, United States**

*Sponsored by*

SPIE

*Cosponsored by*

Becker & Hickl GmbH (Germany) • Carl Zeiss (United States) • Chroma Technology Corp. (United States) • Coherent Inc. (United States) • Semrock Inc. (United States) • ISS, Inc. (United States) • JenLab GmbH (Germany) • Leica Microsystems (United States) • MKS Instruments (United States) • Applied Scientific Instrumentation (United States) • PicoQuant Photonics (United States)

*Published by*

SPIE

**Volume 10498**

Proceedings of SPIE, 1605-7422, V. 10498

SPIE is an international society advancing an interdisciplinary approach to the science and application of light.

The papers in this volume were part of the technical conference cited on the cover and title page. Papers were selected and subject to review by the editors and conference program committee. Some conference presentations may not be available for publication. Additional papers and presentation recordings may be available online in the SPIE Digital Library at [SPIDigitalLibrary.org](http://SPIDigitalLibrary.org).

The papers reflect the work and thoughts of the authors and are published herein as submitted. The publisher is not responsible for the validity of the information or for any outcomes resulting from reliance thereon.

Please use the following format to cite material from these proceedings:

Author(s), "Title of Paper," in *Multiphoton Microscopy in the Biomedical Sciences XVIII*, edited by Ammasi Periasamy, Peter T. C. So, Karsten König, Xiaoliang S. Xie, Proceedings of SPIE Vol. 10498 (SPIE, Bellingham, WA, 2018) Seven-digit Article CID Number.

ISSN: 1605-7422  
ISSN: 2410-9045 (electronic)

ISBN: 9781510614819  
ISBN: 9781510614826 (electronic)

Published by

**SPIE**

P.O. Box 10, Bellingham, Washington 98227-0010 USA  
Telephone +1 360 676 3290 (Pacific Time) · Fax +1 360 647 1445  
[SPIE.org](http://SPIE.org)

Copyright © 2018, Society of Photo-Optical Instrumentation Engineers.

Copying of material in this book for internal or personal use, or for the internal or personal use of specific clients, beyond the fair use provisions granted by the U.S. Copyright Law is authorized by SPIE subject to payment of copying fees. The Transactional Reporting Service base fee for this volume is \$18.00 per article (or portion thereof), which should be paid directly to the Copyright Clearance Center (CCC), 222 Rosewood Drive, Danvers, MA 01923. Payment may also be made electronically through CCC Online at [copyright.com](http://copyright.com). Other copying for republication, resale, advertising or promotion, or any form of systematic or multiple reproduction of any material in this book is prohibited except with permission in writing from the publisher. The CCC fee code is 1605-7422/18/\$18.00.

Printed in the United States of America

Publication of record for individual papers is online in the SPIE Digital Library.

**SPIE. DIGITAL  
LIBRARY**

[SPIDigitalLibrary.org](http://SPIDigitalLibrary.org)

---

**Paper Numbering:** *Proceedings of SPIE* follow an e-First publication model. A unique citation identifier (CID) number is assigned to each article at the time of publication. Utilization of CIDs allows articles to be fully citable as soon as they are published online, and connects the same identifier to all online and print versions of the publication. SPIE uses a seven-digit CID article numbering system structured as follows:

- The first five digits correspond to the SPIE volume number.
- The last two digits indicate publication order within the volume using a Base 36 numbering system employing both numerals and letters. These two-number sets start with 00, 01, 02, 03, 04, 05, 06, 07, 08, 09, 0A, 0B ... 0Z, followed by 10-1Z, 20-2Z, etc. The CID Number appears on each page of the manuscript.

# Contents

- vii *Authors*
- ix *Conference Committee*
- xiii *Introduction*

---

## KEYNOTE SESSION

---

- 10498 04 **Metabolic imaging of tumor for diagnosis and response for therapy (Keynote Paper)**  
[10498-3]

---

## METABOLISM/NADH/FAD/TRYPHOPHAN I

---

- 10498 06 **Ultra-fast HPM detectors improve NAD(P)H FLIM (Invited Paper)** [10498-5]
- 10498 07 **Multiparametric analysis of cisplatin-induced changes in cancer cells using FLIM (Invited Paper)** [10498-6]

---

## METABOLISM/NADH/FAD/TRYPHOPHAN II

---

- 10498 0A **Two-photon luminescence lifetime imaging microscopy (LIM) to follow up cell metabolism and oxygen consumption during theranostic applications (Invited Paper)** [10498-9]
- 10498 0B **Metabolic imaging for breast cancer detection and treatment: a role for mitochondrial Complex I function (Invited Paper)** [10498-10]
- 10498 0G **Multimodal autofluorescence detection of cancer: from single cells to living organism**  
[10498-15]

---

## TECHNOLOGY AND IN VIVO IMAGING I

---

- 10498 0K **Module for multiphoton high-resolution hyperspectral imaging and spectroscopy**  
[10498-19]

---

## FLIM/FRET/FCS I

---

- 10498 0Q **Measuring upconversion nanoparticles photoluminescence lifetime with FastFLIM and phasor plots (Invited Paper)** [10498-25]

---

**FLIM/FRET/FCS II**

---

- 10498 OT **Towards monitoring conformational changes of the GPCR neurotensin receptor 1 by single-molecule FRET (Invited Paper)** [10498-28]
- 10498 OW **FastFLIM, the all-in-one engine for measuring photoluminescence lifetime of 100 picoseconds to 100 milliseconds** [10498-31]
- 10498 OX **In vivo measurement of astrocytic endfoot Ca<sup>2+</sup> and parenchymal vessel responses during 4-AP induced epilepsy using two-photon fluorescence lifetime microscopy** [10498-32]

---

**SECOND/THIRD HARMONIC GENERATION**

---

- 10498 IO **Second harmonic generation microscopy of the living human cornea** [10498-35]

---

**TECHNOLOGY AND IN VIVO IMAGING II**

---

- 10498 I4 **Rapid in vivo vertical tissue sectioning by multiphoton tomography** [10498-39]
- 10498 I6 **In vivo three-photon imaging of deep cerebellum** [10498-41]

---

**COHERENT RAMAN II**

---

- 10498 II **Integrated SRS and fluorescence imaging for study of thermogenesis and lipid metabolism in vivo (Invited Paper)** [10498-52]

---

**COHERENT RAMAN III**

---

- 10498 IP **CARS molecular fingerprinting using a sub-nanosecond supercontinuum light source (Invited Paper)** [10498-59]
- 10498 IQ **In vivo study of lipid synthesis and lipolysis dynamics by stimulated Raman scattering microscopy** [10498-60]

---

**TECHNOLOGY AND IN VIVO IMAGING III**

---

- 10498 IV **In vivo, two-color multiphoton microscopy using a femtosecond diamond Raman laser** [10498-65]
- 10498 IX **Adaptive optics plug-and-play setup for high-resolution microscopes with multi-actuator adaptive lens** [10498-67]

10498 1Y **3D brain oxygenation measurements in awake hypertensive mice using two photon phosphorescence lifetime imaging** [10498-68]

**POSTER SESSION**

---

10498 20 **Simple fibre based dispersion management for two-photon excited fluorescence imaging through an endoscope** [10498-70]

10498 22 **Rapid and lossless bandwidth-switching of a fiber-based optical parametric oscillator for multimodal nonlinear microscopy** [10498-72]

10498 23 **Electronically tunable femtosecond all-fiber optical parametric oscillator for multi-photon microscopy** [10498-73]

10498 24 **Optimizing Ti:Sapphire laser for quantitative biomedical imaging** [10498-74]

10498 25 **Spatially confined photoinactivation of bacteria: towards novel tools for detailed mechanistic studies** [10498-75]

10498 27 **In vivo multiphoton and fluorescence lifetime imaging microscopy of the healthy and cholestatic liver** [10498-77]

10498 28 **Two-photon activation of endogenous store-operated calcium channels without optogenetics** [10498-78]

10498 2A **Comparison of excitation wavelengths for in vivo deep imaging of mouse brain (Student Poster Session Competition)** [10498-110]

10498 2C **Stimulated emission and spontaneous loss pump-probe microscopy for background removal** [10498-81]

10498 2D **Rapid volumetric multiphoton imaging with the combination of an ultrasound lens and a resonant mirror** [10498-82]

10498 2G **Enhancement of measurement speed and photon economy in multiphoton detected fluorescence lifetime imaging microscopy** [10498-85]

10498 2J **Label-free imaging of acanthamoeba using multimodal nonlinear optical microscopy** [10498-87]

10498 2K **Wide field video-rate two-photon imaging by using spinning disk beam scanner** [10498-88]

10498 2L **Improvement of two-photon microscopic imaging in deep regions of living mouse brains by utilizing a light source based on an electrically controllable gain-switched laser diode** [10498-89]

10498 2M **Two-dimensional auto-correlation analysis and Fourier-transform analysis of second-harmonic-generation image for quantitative analysis of collagen fiber in human facial skin** [10498-90]

- 10498 2O **Large scale serial two-photon microscopy to investigate local vascular changes in whole rodent brain models of Alzheimer's disease** [10498-92]
- 10498 2U **Improving multiphoton STED nanoscopy with separation of photons by Lifetime Tuning (SPLIT)** [10498-99]
- 10498 2W **Quantitative 3-dimensional imaging of auxin and cytokinin levels in transgenic soybean and medicago truncatula roots via two-photon induced fluorescence imaging (Student Poster Session Competition)** [10498-101]
- 10498 2Y **Polymer dots enable deep in vivo multiphoton fluorescence imaging of cerebrovascular architecture** [10498-103]
- 10498 2Z **Characterization of a reflective objective with multiphoton microscopy** [10498-104]
- 10498 30 **Improved reference standards for femtosecond three-photon excitation of fluorescence in the wavelength range 950 - 1750 nm** [10498-105]
- 10498 33 **Resolution enhancement of 2-photon microscopy using high-refractive index microspheres** [10498-108]
- 10498 34 **Mueller tensor approach for nonlinear optics in turbid media (Student Poster Session Competition)** [10498-109]
- 10498 36 **Monitoring agrochemical diffusion through cuticle wax with coherent Raman scattering** [10498-113]
- 10498 37 **Novel snapshot hyperspectral imager for fluorescence imaging** [10498-114]